

proaches very closely the "ideal" computer interface with which the users can wire for each channel a single, parallel connector between his instrument(s) and this system and operate his instrument under control of the device linked serially with the STX 1003. Besides operating this interface with teletypewriter and IBM SYS/7, we have also operated it with an Altair 8800 microcomputer. We have used the interface both full duplex and half duplex. We also have a limited capability for selecting the baud rate through control logic.

### Future Serial Computer Interfaces

But despite the versatility of the system built around the STX 1003, microprocessors will undoubtedly play an important role in future serial computer interfaces. With an 8-bit microprocessor and its clock, a parallel I/O port, and a serial I/O port, one could eliminate the need for almost any electronic design. Under program control, one could select and/or change baud rates, codes, numbers of channels transmitted, and formats. Designing microprocessor-controlled serial interfaces would be reduced to wiring a parallel input connector and to programming a PROM (Programmable Read Only Memory). Changing the interface tailored for one instrument to that for another would require changing PROM's. Standard, high-volume serial interfaces between commercially available instruments and commercially available computers would require the same microprocessor-based interface and a much cheaper, mass-produced ROM (Read Only Memory).

In summary, serial computer interfaces can be made quite generally applicable for interfacing many instruments having parallel I/O with computers and many hard-copy devices. Serial interfaces suffer only when required to handle extremely high data rates. Finally, one should not underestimate the power, simplicity, and convenience of a single pair of wires for a computer interface using serial data exchange techniques.

### References

- (1) R. E. Dessy and J. Titus, *Anal. Chem.*, **46**, 294A (1974).
- (2) J. Finkel, "Computer-Aided Experimentation: Interfacing Minicomputers", Wiley, New York, N.Y., 1975.
- (3) D. G. Larsen and P. R. Rony, "The Bugbook IIA", E & L Instruments, Inc., Derby, Conn., 1975.
- (4) D. G. Larsen and P. R. Rony, *Am. Lab.*, **7**, 116 (1975).
- (5) "SERDEX User's Guide", Analog Devices, P.O. Box 280, Norwood, Mass. 02062.
- (6) D. H. Lorimer and A. T. Bell, *Ind. Eng. Chem. Fundam.*, **15**, 71 (1976).

ASK FOR  
YOUR COPY  
TODAY.

**WILKS**

WILKS INFRARED ACCESSORIES

## You'll find 524 IR accessory bargains in this new 40 page catalog.

The 1977 Wilks catalog features accessories for every make and type spectrophotometer. It includes transmission, micro sampling and reflectance accessories, GC/IR and pyrolysis equipment plus accessory kits and crystals. Everything you need for your present or new instrument systems. Quality consistent with Perkin-Elmer, Beckman, and others. You'll save significantly, and get same day order processing from inventory on most items. Call or write today for your copy. Wilks, P.O. Box 449, S. Norwalk, CT. 06856.  
Telephone: (203) 853-1616

**WILKS®**  
INFRARED  
ACCESSORIES

See us at the FACSS MEETING, Booths D-1-4.  
CIRCLE 230 ON READER SERVICE CARD

## The CRYSTAL-CLEAR DROP ...symbol of reliable ULTRA-HIGH-PURITY SOLVENTS

from  
**BURDICK & JACKSON  
LABORATORIES**

Ask for Data Sheet BJ-25  
"Distilled-in-Glass" Solvents  
BJ-21 Liquid Scintillation  
Counting



**BURDICK  
& JACKSON  
LABORATORIES, INC.**

MUSKEGON, MICHIGAN 49442

(616) 726-3171

CIRCLE 27 ON READER SERVICE CARD